




Science and  
Technology  
Facilities Council

# Basic instrument design

Paul F. Henry

Instrument Scientist, ISIS Neutron & Muon Source



# Instrument design steps



## Define a science case!!

Instruments are designed and built to perform defined science

The science case generates the primary requirements


McStas is a tool. If used incorrectly it can be dangerous

# Primary requirements



Science and  
Technology  
Facilities Council

# Defining primary requirements



But also:


- Single crystal, powder, non-crystalline
- Unit cell volume
- Sample environment restrictions
- Need for *in situ* capability
- Sample size / geometry
- Detector type / coverage / stability
- Parasitic scattering
- Etc...

Which come from:


- Science case requirements
- Available technology
- Space limitations
- Budget
- Etc...

# Diffractometer types

$$\lambda = 2d\sin\theta$$



$$\lambda = \frac{2d_c \sin \theta_B}{n}$$



$$\lambda = \frac{3956}{v} = \frac{3956 (t-t_0)}{L_1+L_2}$$


# Diffractometer type

$$\lambda = 2d\sin\theta$$


- Monochromatic (CW)
  - Fix wavelength and scan detector angle
  - Multiple  $2\theta$  required to cover  $Q(d)$  spacing range
  - $Q(d)$  spacing limit  $4\pi/\lambda$  ( $2\pi/d$ )
  - Instrumental count rate factors: Source power, monochromator reflectivity, detector coverage and efficiency, etc
- TOF
  - Fix detector angle and scan wavelength
  - Single  $2\theta$  covers range of  $Q(d)$  space
  - $Q(d)$  range determined by  $\lambda_{\max}$ ,  $\lambda_{\min}$  and  $\theta$
  - Instrumental count rate factors: Source power, moderator performance, beam transport efficiency, detector coverage and efficiency, etc



# Number of reflections



# Reflection density (for CW)





# Resolution

## Monochromatic

$$\frac{\Delta d}{d} = \frac{1}{2} \sqrt{U \cdot \cot^2(\theta) + V \cdot \cot(\theta) + W}$$


- U, V and W are functions of the collimation and U, V also takeoff angle to the monochromator
- Resolution minimum found near the takeoff angle of the monochromator  $2\theta_B$
- Higher takeoff angle gives higher resolution for identical wavelength
- Wavelength produced by monochromator is takeoff angle dependent for any particular hkl plane

## Time-of-flight


$$\frac{\Delta d}{d} = \left[ \Delta\theta^2 \cot^2 \theta + \left( \frac{\Delta t}{t} \right)^2 + \left( \frac{\Delta L}{L} \right)^2 \right]^{1/2}$$

- $\Delta\theta$  is the angular uncertainty
- The main component of  $\Delta t$  is the moderation time of the neutron
- $\Delta L$  is the flight path uncertainty of the neutron mainly due to the finite width of the moderator
- First term can be minimised by moving to higher scattering angle
- Second and third terms minimised by increasing instrument length


# Resolution functions: CW v TOF examples










# Resolution function example: CW



# Resolution function with $\theta_B$ : CW




# Resolution function example TOF: Polaris




Bank 1 – cyan  
Bank 2 – green  
Bank 3 – pink  
Bank 4 – blue  
Bank 5 – orange  
Bank 6 – purple

# Importance of Q-max and Q-resolution



(reproduced from Proffen *et al.* (2003). *Z. Kristallogr.* **218**, 132-143) showing Ni  $G(r)$  from 2 instruments with different Q-resolution with  $Q_{\text{max}}$  cut-off of  $35 \text{ \AA}^{-1}$




- Good  $\Delta d/d$  resolution at low  $d$  ( $\Delta Q/Q$  resolution at high  $Q$ )
- Access to as low  $d$  (high  $Q$ ) as possible

# Summary of Q resolution and range

- CW:
  - Simple, symmetric peakshape function
  - Best resolution where diffraction peak density is highest in scattering angle
  - Different wavelength can be used to give Q resolution where required
  - Different takeoff angle can be used to change resolution function and wavelength
  - Instrument can have high Q resolution over a very limited Q range
- TOF:
  - Complex asymmetric peakshape related to moderator characteristics
  - Instrument length and moderator give wavelength band and overall resolution
  - Q resolution almost constant for a given detector bank so increasing peak density with Q can be an issue
  - Q resolution improved by moving to higher scattering angle detector bank
  - Q range determined by scattering angle of detector bank



# Defining primary requirements



But also:

- Single crystal, powder, non-crystalline
- Unit cell volume
- Sample environment restrictions
- Need for *in situ* capability
- Sample size / geometry
- Detector type / coverage / stability
- Parasitic scattering
- Etc...

Rank them in importance to the proposed instrument role




# Instrument layout



Science and  
Technology  
Facilities Council

# Instrument layout: follow the neutrons



Understand the impact of complexity on performance v cost

Identify potential issues early

Identify necessary development / proof of concept studies


Prepare yourself for discussions with engineers and vendors

# Source



Science and  
Technology  
Facilities Council

# Source influence on CW or TOF




# Moderators


	cold	thermal	hot
moderator	liquid D <sub>2</sub>	Liquid D <sub>2</sub> O	graphite
moderator temperature	20K	300K	2000K
neutron wavelength	3→20Å	1→3Å	0.3→1Å
sample lengthscale	1Å→100 nm	0.3→5Å	0.1→2Å
sample timescale	1kHz→1 THz	0.1→10 THz	1→100 THz

- Wavelength range
- Peak-shape (residence time)
- Temperature stability
- Viewing angle
- Viewing multiple moderators
- Poisoning
- Moderator development


# Moderators: short pulse



# Moderators: long pulse






# Source brightness for CW or TOF





# Source stability / availability



Spallation sources have lower availability than reactors. Does this compromise science case?

# Instrument type v source type

- Reactors build CW instruments\*
  - Low peak brilliance, high time-average brilliance
  - Variable reflectivity from monochromators limit low  $\lambda$  use
  - High Q not easily reached
  - Match moderator and monochromator take-off angle to Q range and resolution
  - Beam always on

\*Except when significantly restricted geometry constraints from science case necessitate use of TOF

- Pulsed sources build TOF instruments#
  - High peak brilliance, low time-averaged brilliance
  - Require efficient beam transport
  - High Q possible
  - Increase instrument length to improve resolution at expense of bandwidth
  - Variable Q range and resolution from detector angles
  - Beam availability can compromise science


#Remains to be seen for long pulse sources



Science and  
Technology  
Facilities Council

# Summary

Q(d) Resolution



Q(d) range

Count-rate

- Different source types offer complementary possibilities
- Don't limit yourself to a particular source type without good reason



It may be that our spallation source is not ideally suited to your science case requirements!

# Beam transport




Science and  
Technology  
Facilities Council

# Monochromator reflectivity for CW




Choose a material and plane to access necessary Q-range  
Reflectivity falls as wavelength decreases

# Q-range accessible with CW




# Higher reflection order contamination: CW






# Beam transport for pulsed sources

- Use distance to separate wavelengths
- Need choppers to prevent frame overlap
- Moderator pulse-width and distance determine resolution
- Source repetition rate determines available time-window
- Beam transport as a function of wavelength crucial











# Q-range with TOF (available bandwidth)




# Q-range with TOF (available bandwidth)



# Q-range with TOF (available bandwidth)



# Q-range with TOF (available bandwidth)



# Q-range and resolution for CW and TOF

For CW:

- For monochromatic instruments the  $Q_{\max}$  is  $4\pi/\lambda$  i.e. when  $\sin\theta = 1$ ,  $\theta = 90^\circ$ ,  $2\theta = 180^\circ$
- If a high  $Q_{\max}$  is required a shorter wavelength must be used.
- Shorter wavelengths are produced by higher order hkl planes
- Reflectivity is lower for shorter wavelengths
- Realistic  $Q_{\max}$  of around  $25 \text{ \AA}^{-1}$

For TOF:

- $Q_{\max}$  depends on  $\lambda_{\min}$  and detector  $\theta$ .
- $\lambda_{\min}$  can be much lower than for the CW case allowing  $Q_{\max} > 100 \text{ \AA}^{-1}$
- $\lambda_{\min}$  determined by the moderator, transport characteristics of the guide and which frame the instrument is working in



# Other factors

- Line-of-sight loss v  $T_{\text{zero}}$  v beamstop shielding v background v etc...
- Phase space homogeneity (primarily divergence)
- Footprint of guide
- Need for /positioning of choppers
- Shielding requirements
- £ – longer = more expensive
- £ – higher M-coatings = more expensive
- £ – complex guide geometries = more expensive
- £ – large guide cross-sections = more expensive
- £££!



Sample area



Science and  
Technology  
Facilities Council





# Vacuum v air?









# Restricted geometry: Pearl



Instrument geometry designed to match sample environments

# Parasitic scattering: collimation

NaAlCaF 5mm, 2.4AA HR (10' 5/5), TasMag ROC



# Other factors


- Beam definition close to sample – slits and / or guide snouts
- Accessing sample environment during experiment
- Sample environment scattering / materials
- Difficulty running under vacuum v background / beam attenuation from air scattering
- Shadowing from collimation
- Sample geometry
- Non-magnetic components to allow magnetism studies
- Future upgrades – e.g. polarisation
- Etc...

# Detectors




Science and  
Technology  
Facilities Council

# Detector efficiency



$^3\text{He}$  detection efficiency as a function of detection depth



Predicted detector efficiency CASCADE-detector for 20  $^{10}\text{B}$  layers

# Detector coverage v sample geometry

Appl Cryst

JAC

Journal of Applied Crystallography

IUCr

IT

WDC

search IUCr Journals

[home](#) [archive](#) [editors](#) [for authors](#) [for readers](#) [submit](#) [subscribe](#) [open access](#)

JAC

RESEARCH PAPERS

*J. Appl. Cryst.* (2011). **44**, 295-298  
<https://doi.org/10.1107/S0021889811005498>  
Cited by [2](#)



## The ABC of powder diffractometer detector coverage

**K. H. Andersen, P. M. Bentley and L. D. Cussen**

This article addresses the question of the most effective detector configuration for neutron spectrometers using a discussion of a particular case - constant-wavelength powder diffractometers at continuous neutron sources. A first variation uses an essentially one-dimensional 'banana' detector coupled with out-of-scattering-plane beam divergence before and after the sample. A second variation uses an incident beam tightly defined both in- and out-of-plane coupled to a '4 $\pi$ ' detector covering all possible scattering angles after the sample. It is widely believed that the 4 $\pi$  arrangement is superior for most varieties of neutron spectrometer but is more difficult and more expensive to implement. Starting from a commonly used overall instrument figure-of-merit, this article presents simple arguments leading to the surprising conclusion that this is untrue for these constant-wavelength powder diffractometers, provided only that the cylindrical samples used with the banana detector have a height greater than 2.4 times their diameter.

**Keywords:** neutron powder diffraction; constant-wavelength powder diffractometers; detector configuration.

[Read article](#) [Similar articles](#)

ACCESS



Follow J. Appl. Cryst.

E-alerts

Twitter

Facebook

RSS

Search IUCr Journals

doi

Author

All journals

volume

page

Advanced search

Copyright © International Union of Crystallography  
Home Contact us Site index About us Partners and site credits Help Terms of use Privacy and cookies  
The IUCr is a scientific union serving the interests of crystallographers and other scientists employing crystallographic methods.



Science and  
Technology  
Facilities Council

# Other factors

- Count-rate / saturation
- Known / proven technology?
- Pixel size
- Sample geometry / shadowing
- Performance variation
- Static or moveable?
- Restricted geometry
- Detector stability v time v T
- Maintenance / replacement
- In vacuum?
- Unit £





# Figure of Merit



Science and  
Technology  
Facilities Council



# What to use in the figure of merit?

- Flux at sample position
  - Total count-rate on detector
  - Divergence profile
  - Brilliance transfer
  - Resolution of a particular diffraction peak(s) at detector
  - Number of unwanted neutrons
  - Combinations
- 
- Effects of sample volume / geometry
  - Detector layout / varying resolution



# McStas optimisation



Science and  
Technology  
Facilities Council

# What might be optimised in McStas?

- Detector layout in plane
- Detector coverage out of plane
- Sample size / geometry
- M-coatings of a complex guide
- Cost v performance of a guide
- Divergence acceptance / transport
- Chopper positioning / performance
- Maximising good:bad neutron transport
- Optimal detector resolution / pixel size vs sample geometry vs beam transport
- Etc...









Science and  
Technology  
Facilities Council

# Final thoughts

# Instrument design v reality: GEM






# Instrument design: Polaris rebuild v 2006






# Instrument design: Polaris










# Instrument reality: Polaris





# Instrument performance: Polaris



Bank 1 – cyan  
Bank 2 – green  
Bank 3 – pink  
Bank 4 – blue  
Bank 5 – orange  
Bank 6 – purple

- Increased count rate  $\times 3$  at high scattering angle to  $>20$  for low angle banks
- Resolution improvement e.g. bank 5 and 6 of  $3 \times 10^{-3}$  cf.  $5 \times 10^{-3}$
- Improvement in data at high Q



Science and  
Technology  
Facilities Council

The background features a large blue rectangle on the right side, which is partially overlaid by a series of blue lines forming a complex, abstract geometric pattern. This pattern includes several nested, right-angled shapes and lines that create a sense of depth and movement. The overall color scheme is dominated by orange and blue.

# Questions?