


Erik Knudsen
Peter Willendrup
Emmanuel Farhi

And many other
contributors:

Manuel Sanchez del Rio
Claudio Ferrero
Kell Mortensen
Robert Feidenshans'l
Kim Lefmann

McXtrace



McXtrace: X-ray Monte-Carlo ray-tracing



What is McXtrace ?

A software to simulate X-ray beam lines with incoherent and coherent sources, including optics, samples (scattering), detectors.

This is a ray-tracing code (Monte-Carlo).

Can easily couple to other codes.

Made to be extensible. Massively parallel.



"McXtrace: a Monte Carlo software package for simulating X-ray optics, beamlines and experiments",
Journal of Applied Crystallography, vol. **46** (2013) 679.





What can it be used for ?

Understand beam line as a whole.

Improve beam lines. “What if the ring lattice is changed ?”


Train users – plan experiments.

Help analyse data (limited by sample models).



How does it work

1. Describe a beam-line in the McXtrace language (text file).
2. Automatically convert beamline into ANSI C
3. Compile (Automatically)
4. Run (Automatically)




1. Optimized for your platform
2. Only includes what you use

Built over proven McStas (neutrons) technology (1998-).

Produce data sets and plots, as well as 3D views.

Typical simulation time: 10s, up to few minutes.

User interfaces are in Python (also keep legacy Perl). Matlab interface.






A photon particle is...

Photon ray/package:

- $(\mathbf{r}, \mathbf{k}, \phi, t, p, \mathbf{E})$
- \mathbf{r} - spatial coordinates
- \mathbf{k} - wave vector
- ϕ - phase
- t - time
- p - photon weight
- \mathbf{E} - Electrical field polarisation





Installation

Available on all systems (Windows, Linux, Mac) at <<http://www.mcxtrace.org/>>.

[1.5 RC1: <http://downloads.mcxtrace.org/mcxtrace-1.5rc1/>]

Available as Virtual Machine <<http://ord03246.synchrotron-soleil.fr>> (exp.)

Available example at <<https://sim.e-neutrons.org>>



Help / Community

Extensive manuals, tutorials, examples.

Forum/mailing list.

Built over McStas – vast knowledge, community and expertise.

Code hosted at <<https://github.com/McStasMcXtrace/McCode>>





What you get in the box



A software that simply runs...

- **Components (80)**
 - Sources
 - Optics
 - Samples
 - Detectors
- **A GUI to assemble and run BL.**
- **Tools to view the results (e.g. detectors).**
- **Tools to view the BL geometry.**
- **Examples (40)**
- **Extensive documentation, mailing list...**



User interfaces *and* scripting (bash, Python, Matlab).

Single shot and parameter-scanned simulations.

Generated data files as text or HDF5/NeXus.

Fully adapted to HPC (MPI).

Optimisation (swarm, GA, simplex, Hooke, Powell, ...) via *iFit* <<http://ifit.mccode.org>>





Sources

Synchrotron ID


- Bending magnet [B.D. Patterson, *Am. J. Phys.* **79**, 1046 (2011)]
- Undulator [K.J. Kim, *AIP, conf. proc.*, 184, 1989]
- Wiggler [B.D. Patterson, *Am. J. Phys.* **79**, 1046 (2011)]

Lab/ideal stuff

- Laboratory X-ray tube (e.g. rotating anode)
- Ideal, point and gaussian

Interfaces with other software

- **Spectra** (R) <<http://spectrax.org/spectra/>>
- **Simplex** (R) <<http://spectrax.org/simplex/index.html>>
- **Genesis** (R) <<http://genesis.web.psi.ch/>>
- **Shadow** (RW) <<https://github.com/oasys-kit/shadow3>>
- **MCPL** (GEANT4, PHITS, MCNP) (RW) <<https://mctools.github.io/mcpl/>>






Optics


Enough to start with:

- Bragg crystal (monochromator, incl. bent)
- Capillary
- Filter (absorption and refraction)
- Lenses
- Mirrors (flat, curved, multi-layers, twin KB multi-layer)
- Zone plate
- Grating (lamellar, blazed)
- Slit, beam-stop, ...



More to come (with your help ?)



Components can be arranged in groups.





Samples

- SAXS sample (60 models from SasView, PDB, Nanodiscs, Liposomes, $I(q)$, ...)
- Powder (diffraction)
- Polycrystal (diffraction)
- Pump-probe (2 states) molecule
- Single crystal (diffraction – **hold on !**)



All samples can have simple geometric shapes (incl. hollow).
Powder and SX can have any shape (PLY/STL).
Powder sample supports multiple concentric geometries (e.g. for cryostat).
McXtrace comes with a material data base, and can use e.g. NIST files.





Detectors

Currently only record distributions (1D, 2D, event lists, ...) for any combination of state parameters (position, divergence, energy, power, phase, E-field, ...).

Not satisfactory !



Actual detector efficiency to be added (simple models exist), e.g.

- <<https://bl831.als.lbl.gov/~jamesh/mlfsom/>>



Examples

SAXS

Powder

Single crystal

Coherence

MaxII 711 XRD/MX

DanMAX@MaxIV Powder XRD 2D

https://sim.e-neutrons.org/instrument-menu/xrd/MAXIV_DanMAX_pxrd2d



SAXS

```

1 /*****
2 *      McXtrace instrument definition URL=http://www.mcxtrace.org
3 *
4 * Instrument: TestSAXS
5 *****/
6
7 DEFINE INSTRUMENT TestSAXS(
8   DistanceFromSourceToFirstPinhole = 0.05, DistanceFromSourceToSecond
   DistanceFromSecondPinholeToSample = 0.6, DistanceFromSampleToDetector
   RadiusOfDetector = 0.1, Lambda = 1.54, DLambda = 0.01, SAMPLE=0 )
9
10 TRACE
11 COMPONENT Source = Source_flat(
12   xwidth = 0.0005, yheight = 0.0005, dist = DistanceFromSourceToSecond
   focus_yh = 0.0003, focus_xw = 0.0003, lambda0 = Lambda, dlambd = DL
13 AT (0, 0, 0) ABSOLUTE
14
15 COMPONENT FirstPinhole = Slit(radius = 0.0002)
16 AT (0, 0, DistanceFromSourceToFirstPinhole) RELATIVE Source
17
18 COMPONENT SecondPinhole = Slit(radius = 0.00015)
19 AT (0, 0, DistanceFromSourceToSecondPinhole) RELATIVE Source
20
21 COMPONENT Sample0 = SAXSSpheres(
22   xwidth = 0.01, yheight = 0.01, zdepth = 0.01, R = 50.0, SampleToDetectorDistance =
   DistanceFromSampleToDetector, DetectorRadius = RadiusOfDetector )
23 AT(0, 0, DistanceFromSampleToDetector) RELATIVE SecondPinhole
24
25 COMPONENT Beamstop = Beamstop(
26   radius = 0.001)
27 AT (0, 0, DistanceFromSampleToDetector - 0.000001) RELATIVE Sample0
28
29 COMPONENT Monitor = PSD_monitor(
30   filename = "PSDMonitor", xwidth = 2.0 * RadiusOfDetector / sqrt(2
   * RadiusOfDetector / sqrt(2.0), nx = 200, ny = 200, restore_xray =
31 AT (0, 0, 0.000001) RELATIVE Beamstop
32
33 COMPONENT QMonitor = SAXSQMonitor(
34   RadiusDetector = RadiusOfDetector, DistanceFromSample =
   DistanceFromSampleToDetector, LambdaMin = Lambda, Lambda0 = Lambda,
   2000 )
35 AT (0, 0, 0.0000010000000001) RELATIVE Beamstop
36
37 END
38

```

Start simulation

Instrument parameters (D=floating point, I=integer, S=string)

DistanceFromSourceToFirstPinhole: 0.05	DistanceFromSourceToSecondPinhole: 0.7	DistanceFromSecondPinholeToSample: 0.6
DistanceFromSampleToDetector: 0.48	RadiusOfDetector: 0.1	Lambda: 1.54
DLambda: 0.01	SAMPLE: 0	

Simulation

Simulation/Trace: Simulation

Autoplot: - None --

Ray count: 1000000

Output subdir (optional):

Sweep steps (optional):

MPI: No clustering

MPI node count: 2


Advanced

Random seed:

Gravity: Off

Start Cancel

11 s



Powder Diff.

```

1 /*****
2 * Instrument: Test_PowderN
3 * %Parameters
4 * TTH: [deg.] Two theta rotation for detector arm
5 * %End
6 *****/
7
8 DEFINE INSTRUMENT Test_PowderN(TTH=13.4)
9
10 TRACE
11 COMPONENT src = Source_flat(
12     yheight = 1e-3, xwidth = 1e-3, dist = 10, focus_xw = 1e-3,
13     focus_yh = 1e-3, E0 = 45, dE = 1)
14 AT (0, 0, 0) ABSOLUTE
15
16 COMPONENT psd0 = PSD_monitor(
17     nx = 100, ny = 100, filename = "psd0", xwidth = 2e-3, yheight = 2e-3)
18 AT (0, 0, 1e-9) RELATIVE src
19
20 COMPONENT sample = PowderN(
21     reflections = "Fe_bcc229_lt13deg.dat", format=Crystallographica,
22     material = "Fe.txt", radius = .5e-4, yheight = 1e-3, pack = 0.5, Vc =
23     p_inc = 0, p_transmit = 0.1, DW = 0, d_phi = 45)
24 AT (0, 0, 10) RELATIVE PREVIOUS
25 EXTEND %{
26     if(!SCATTERED) ABSORB;
27 }%
28
29 COMPONENT psd4pi = PSD_monitor_4PI(
30     nx = 180, ny = 180, filename = "psd4pi", radius = 0.1, re
31 AT (0, 0, 0) RELATIVE sample
32
33 COMPONENT detector = PSD_monitor(
34     nx=200, ny=200, xwidth=0.1, yheight=0.1, filename="psd1",
35 AT(0,0,0.1) RELATIVE sample
36
37 COMPONENT ttharm= Arm()
38 AT(0,0,0) RELATIVE sample
39 ROTATED (0,TTH,0) RELATIVE sample
40
41 COMPONENT detector2 = PSD_monitor(
42     nx=200, ny=200, xwidth=2e-3, yheight=1e-2, filename="psd2",restore_xray=1)
43 AT(0,0,0.1) RELATIVE ttharm
44 END

```

Start simulation

Instrument parameters (D=floating point, I=integer, S=string)

TTH:

Simulation

Simulation/Trace:

Autoplot:

Ray count:

Output subdir (optional):

Sweep steps (optional):

MPI:


MPI node count:

Advanced

Random seed:

Gravity:

6 s



Single crystal diff.

Start simulation

Instrument parameters (D=floating point, I=integer, S=string)

TTH:

Simulation

Simulation/Trace:

Autoplot:

Ray count:

Output subdir (optional):

Sweep steps (optional):

MPI:



MPI node count:

Advanced


Random seed:

Gravity:

↓ 32 s



Very early stage !






2 Slits interferences

```

1 /*****
2 * Instrument: template_2Slit_Diff
3 * An example instrument showing a Young's double slit experiment
4 * %Parameters
5 * SLITW: [m] Width of the slit in question
6 * %End
7 *****/
8
9 DEFINE INSTRUMENT template_2Slit_Diff(SLITW=1e-6, SLITSEP=1e-5)
10
11 DECLARE %{
12     double slit_offset;
13 %}
14
15 INITIALIZE %{
16     slit_offset=SLITSEP*0.5;
17 %}
18
19 TRACE
20
21 COMPONENT Origin = Progress_bar()
22     AT (0,0,0) ABSOLUTE
23
24 COMPONENT source=Source_pt(focus_xw=6e-6, focus_yh=0.8e-6, dist=1,
25     lambda0=2.0, dE=0, gauss=1, randomphase=0)
26 AT(0,0,0) RELATIVE Origin
27
28 COMPONENT s1 = Slit(
29     xwidth=SLITW, yheight=0.8e-6, dist=1, focus_xw=200e-6, focus_yh=2e-6,
30     focus_x0=0.0, focus_y0=0.0 )
31 AT(slit_offset,0,1) RELATIVE source
32 GROUP slits
33 COMPONENT s2 = COPY(s1)()
34 AT(-slit_offset,0,1) RELATIVE source
35 GROUP slits
36
37 COMPONENT psd0 = PSD_monitor_coh(
38     yheight=2e-6, xwidth=200e-6, nx=501, ny=1, filename="psd0")
39 AT(0,0,2) RELATIVE source
40
41 END

```

↓ 1 s






Max II 711 (MX)


A complete BL in 200 lines.

Execution: 40 s

3D view with 'mxdisplay'
(WebGL in a browser)



Sample:
Diamond
powder






Max IV DanMAX (PXRD)

Try it on the web:

https://sim.e-neutrons.org/instrument-menu/xrd/MAXIV_DanMAX_pxrd2d

A complete BL in 500 lines.

Execution: **30 s**





Limitations

Initially oriented towards 'science', not machine physics. That's why we favour interfaces with other packages.

Can handle 100% coherent OR 100% incoherent beam.

BUT

Made to ease contribution and extend. Easy to learn.

Possibility to add new optics, samples and detectors. This is Open Source.

Particle (photon) description can be extended easily (add more state parameters to e.g. better describe coherent beams) to handle e.g. polarised BL.





The Future of McXtrace

Contribute to new sample kernels:

- Single Xtal (validate)
- EXAFS
- IXS/XPS/ARPES [$S(\mathbf{q}, w)$ exists]

Contribute to new detector models:

- From MLF-SOM (= MOS-FLM⁻¹)
- Interfaces with e.g. HORUS, GEANT4, ...

Add more optics.

Tomography, ptychography...

On-going work to use GPU's *via* NVIDIA OpenACC.

Add interface to OpenPMD (HDF5) particle/ray records

Couple to e.g. OASYS, SIMEX, COMSYL, SRW, ... SpotX, Solemio ?

Tutorial/School as a Training for Soleil Staff planned at Soleil: **2-6 December 2019**





How to contribute ?

Using McXtrace is already a good start. You can get a feeling in a few minutes...

You may then build new beam line descriptions. Generally 100-200 lines.

You further may adapt or build new components, e.g. 100-200 lines.

Send us an email or push on GitHub to allow others to benefit from your work:

emmanuel.farhi@synchrotron-soleil.fr
erkn@fysik.dtu.dk

<http://www.mcxtrace.org/>