

Mantid – McStas Exercise

Torben Nielsen
Peter Willendrup



www.europeanspallationsource.se

May 2021

Exercise 1:


<https://github.com/McStasMcXtrace/McCode/wiki/McStas-and-Mantid>

- Reproduce the SANS results shown on the wiki page in section 8: "MantidPlot view of McStas event data" and section 9: "Mantid reduction of McStas event data"
- **Step one:** Convert templateSANS.instr to a Mantid compatible version templateSANS_Mantid.instr
- **Step two:** Run the simulation and load data into Mantid
- **Step three:** Derive $I(q)$ and compare to SasView data



Exercise 2

- Reproduce the SANS results shown on the wiki page in section 10: "ISIS SANS2D"
- Go to the wiki-page:
<https://github.com/McStasMcXtrace/McCode/wiki/McStas-and-Mantid#isis-sans2d>
- Reproduce data shown in Figure 8:



- Copy folder:
 - `/nfs/www/html/users/troland/nexus_local`
 - to your home dir
-
- `cp -r /nfs/www/html/users/troland/nexus_local .`
 - `cd nexus_local`
 - `more how-to-run.txt`

Lots of text in how-to-use.txt file

- # compile for mpi
- run
- # Mantid
- open new terminal
- go to same node
- # Load Mantid
- module load mantid/5.1.1
- # mantid reduction
- mantidpython MantidQ.py
- # Back on McStas client
- python plot_lq.py
- # use model nr5: spheres 25 AA
- # use model nr6: spheres 500 AA
- # use model nr15: spheres 150 AA
- # change pinhole S6,A2 -> 0.001
- # do we see any change