

# Optimisation of NR and SANS Instrumentation




Science & Technology Facilities Council

ISIS



- Length scales probed - 1s to 100s nm
  - Studies include surfactants, polymers, liquid crystals, nanoparticles, lipids and fibres.



## Science and Sample Environment

**Solution scattering**  
Sample Changer

**Exchange in nanoemulsions**  
Stopped-Flow

**Defects in metals**  
Sample Changer

**Foams**  
Foam Cell

**Hydrogen loading**  
Gas Handling

**Flux line lattices**  
17T cryomagnet

**Polymer mixtures**  
T-jump cell

**Skymions**  
3D magnet

**Templating of nanoparticles**  
Sample Changer

**Colloidal crystals**  
Mini Huber Stack

**Interfacial structures of polymers  
at various interfaces**  
GISANS

**Organic Light Emitting Diodes**  
GISANS

**Fuel cell membranes**  
Humidity Chamber

**Drug Delivery Systems Under Shear**  
Rheometer

**Surfactant mixing**  
Microfluidics set up

**Nanoparticles in metal alloys**  
Furnace and 1.5T magnet

**Movement of drugs through  
and into vesicle bilayers**  
Stopped-Flow

**Growth of fibrils**  
Stopped-Flow

**Micellization in CO<sub>2</sub>**  
CO<sub>2</sub> Pressure cell

**Food products**  
Shear cell

**Magnetic nanoparticles**  
1.5T magnet

**Structural colour**  
Sample Changer

**Cryoprotectants**  
Linkam

**Interaction of polymers with DNA**  
Sample Changer


**Contrast variation**  
NURF set up

**(SE)SANS**

```
graph TD; SANS((SE)SANS) --> SolutionScattering[Solution scattering<br/>Sample Changer]; SANS --> Exchange[Exchange in nanoemulsions<br/>Stopped-Flow]; SANS --> Defects[Defects in metals<br/>Sample Changer]; SANS --> Foams[Foams<br/>Foam Cell]; SANS --> Hydrogen[Hydrogen loading<br/>Gas Handling]; SANS --> Flux[Flux line lattices<br/>17T cryomagnet]; SANS --> Polymer[Polymer mixtures<br/>T-jump cell]; SANS --> Skymions[Skymions<br/>3D magnet]; SANS --> Templating[Templating of nanoparticles<br/>Sample Changer]; SANS --> Colloidal[Colloidal crystals<br/>Mini Huber Stack]; SANS --> Interfacial[Interfacial structures of polymers<br/>at various interfaces<br/>GISANS]; SANS --> OLED[Organic Light Emitting Diodes<br/>GISANS]; SANS --> FuelCell[Fuel cell membranes<br/>Humidity Chamber]; SANS --> DrugDelivery[Drug Delivery Systems Under Shear<br/>Rheometer]; SANS --> Surfactant[Surfactant mixing<br/>Microfluidics set up]; SANS --> Nanoparticles[Nanoparticles in metal alloys<br/>Furnace and 1.5T magnet]; SANS --> Movement[Movement of drugs through<br/>and into vesicle bilayers<br/>Stopped-Flow]; SANS --> Fibrils[Growth of fibrils<br/>Stopped-Flow]; SANS --> Micellization[Micellization in CO2<br/>CO2 Pressure cell]; SANS --> Food[Food products<br/>Shear cell]; SANS --> Magnetic[Magnetic nanoparticles<br/>1.5T magnet]; SANS --> Structural[Structural colour<br/>Sample Changer]; SANS --> Cryo[Cryoprotectants<br/>Linkam]; SANS --> Polymers[Interaction of polymers with DNA<br/>Sample Changer]; SANS --> Contrast[Contrast variation<br/>NURF set up];
```



- Lengthscales are explored in reciprocal space by detecting the number of scattered neutrons as a function of the scattering vector,  $Q$ .



Units are either  $\text{\AA}^{-1}$  or  $\text{nm}^{-1}$  i.e. the smaller the value of  $Q$  the bigger the object

$Q$  is also related to wavelength and the scattering angle by:


$$Q = \frac{4\pi \sin\left(\frac{\theta}{2}\right)}{\lambda}$$

$Q$  (size) range is varied by altering  $\theta$  or  $\lambda$

- $L_1 = L_2$  for optimal  $Q$  resolution
- To reach the smallest  $Q$  values the incident flux is always lower in conventional 'pinhole collimation' SANS as a long incident collimation is needed



- At ISIS we use time-of-flight (TOF) to record numerous diffraction patterns at different wavelengths, then we combine the patterns at different wavelengths onto a single Q scale
- TOF SANS therefore has a wide Q range in a single measurement



On Sans2d:


- Wavelength range = 1.75 - 16.5 Å
- 2 x 1 m<sup>2</sup> detectors
- Q-range = 0.001 - 3 Å<sup>-1</sup> - dependent on detector positions

Need wavelength dependent corrections for:


- Monitor spectrum
- Detector efficiency
- Sample transmission (measured)



- Probe thin films at an interface which can be buried under neutron transparent material.
- Analogous effect to the rainbow generated by light reflecting from an oil film on water. Reflection and refracted beams are observed.
- $\theta_i = \theta_f$  for specular reflection which probes the direction perpendicular to the interface




$$Q = \frac{4\pi \sin(\theta)}{\lambda}$$




Incident Beam

Reflected Beam






- Probes both structure perpendicular to the interface and in the plane of the sample along the beam direction.
- McStas doesn't have a sample to model this.






- $Q_z$  Perpendicular to the interface
- $Q_x$  In plane structure parallel to the beam
- $Q_y$  In plane structure perpendicular to the beam





# Neutron reflection science themes



# Instrument Optimisation

For NR and SANS the key to everything is background.

NR signals  $10^{-7}$  or lower

SANS  $< 0.1 \text{ cm}^{-1}$


McStas is not the best tool to determine this.



Science & Technology Facilities Council

ISIS

# Larmor vs. SANS2D



Larmor – Open circles  
SANS2D – Filled circles

Empty instrument  
backgrounds scaled to  
absolute units.

Factor of >3 background  
increase on Larmor. For  
weak signals (<0.1) a  
factor of 10 increase in  
counting time is needed  
to extract data with  
similar uncertainties






# Science Focus

- SANS



- Science focus determines instrument requirements and requires a science case. For example:

- Broad simultaneous  $q$  range
    - Smallest  $q$ s
    - High flux with relaxed resolution
    - Polarisation
    - Complex sample environments
    - Grazing Incidence SANS



# Science Focus

- NR
  - Science focus strongly determines instrument requirements
    - Broad simultaneous  $q$  range
    - Small samples
    - Specular reflection only
    - Polarisation
    - High flux with relaxed resolution
    - Complex sample environments





# Groundwork



- Basic questions to be answer with spreadsheets
  - Q range required
  - Possible Wavelength Band
  - Chopper positioning
    - Need to know what is technically possible
  - Detector options
    - Pixel size and coverage





# Groundwork 2




- Other limitations:  
Larmor
  - Many physical constraints



- Physical constraints strongly limited design options

# Performance Metrics

- Unlikely to be just one
  - Flux on sample – a starting point only
  - Homogeneous phase space transport
  - Samples will be needed






# NR Example

- Minimise the data collection time required to measure a 40x40mm sample consisting of a 15Å thick layer of SiO<sub>2</sub> on a silicon block against a pool of D<sub>2</sub>O.
  - Background is ignored.
  - Q resolution must be considered
  - What q value is it important to reach?
  - What about thicker films?







a) Si-15ÅSiO<sub>2</sub>-air.txt  
 $\Delta\theta/\theta=4\%$ , WFM OFF




b) Si-15ÅSiO<sub>2</sub>-air.txt  
 $\Delta Q/Q=4.5\%$ , WFM ON





SiO<sub>2</sub> on Si  
air

a) air-Si-15ÅSiO<sub>2</sub>-300ÅdPS.txt  
 $\Delta\theta/\theta=4\%$ , WFM OFF



b) air-Si-15ÅSiO<sub>2</sub>-300ÅdPS.txt  
 $\Delta Q/Q=4.5\%$ , WFM ON



dPS on Si  
air



- SANS – There are lots:
  - The 3 most realistic/useful are
    - SANS\_benchmark2
    - SasView\_model
    - SANS\_spheres2
- NR - Just one:
  - Multilayer\_Sample

# SANS\_Benchmark2

COMPONENT sample = SANS\_benchmark2(  
xwidth=0.01,  
yheight=0.01,  
zthick=0.001,  
model=model\_nr, - model number to be picked from the list  
dsdw\_inc=0.0, - phenomenological incoherent background  
sc\_aim=1.0, - fraction to scatter (the rest transmit)  
sans\_aim=1.00, - fraction to SANS scatter the rest are incoherent  
singlesp = 1.0) - Single or multiple scatter  
AT (0, 0, 0.0) RELATIVE sampleMantid




# Multilayer\_Sample

COMPONENT samp1 = Multilayer\_Sample(  
xwidth = 0.05, - width of sample  
zlength = 0.15, - length of sample along the beam  
mu\_inc=0.138, - incoherent scattering length  
ythick=0.01, - thickness of the sample subphase  
nlayer = 0, - number of layers  
sldPar = {0.0,6.35e-6}, - list of layer scattering length densities  
dPar = {0.0}, - list of layer thicknesses  
sigmaPar = {5.0}, - list of interface roughnesses  
target\_index=1, - target to focus the incoherent scattering onto  
focus\_xw=2\*tend, - width of incoherent scattering focus  
focus\_yh=0.01, - height of incoherent scattering focus  
frac\_inc=FRAC) - fraction of neutrons to scatter incoherently  
AT (0.0, spos\*t15, spos) ABSOLUTE



# Chopper Diagrams



A simple spreadsheet will help with chopper positioning and will show where frame overlap problems will occur.

$1\text{\AA}$  neutron velocity =  $3956.0339760560055\text{ ms}^{-1}$

<https://webapps.frm2.tum.de/neutroncalc/>




Science & Technology Facilities Council

ISIS

# $T_0$ Chopper

- Chopper to block the beam at time 0 when a flash of high energy neutrons passes down the beamline.




- How do you model this?



# Q Resolution - SANS

- Mildner & Carpenter J.Appl.Cryst. 17(1984)249-256
- Assuming, small angles.
- Circular apertures, isotropic scattering,  $\Delta R$  (Detector resolution) and  $\Delta\lambda$  are assumed to be rectangular distributions not Gaussians.
- Optimal SANS configuration when  $L_1=L_2$  and  $R_1=2R_2$
- Weak constraint: Changing  $L_1$  and  $L_2$  allows the wavelength band be changed





$$\left(\frac{\sigma_Q}{Q}\right)^2 = \left(\frac{R_1 L_2}{2R L_1}\right)^2 + \left(\frac{R_2 (L_1 + L_2)}{2R L_1}\right)^2 + \frac{1}{12} \left(\frac{\Delta R}{R}\right)^2 + \frac{1}{12} \left(\frac{\Delta \lambda}{\lambda}\right)^2$$





# Phase Space

- Homogeneous illumination of the sample is important for normalisation in both NR and SANS
  - Add a PSD monitor after A2 to probe the beam shape before the sample
  - Add 2 divergence monitors to probe the divergence transported to the sample in both X and Y



# ISIS Moderator Component

- The new component and moderator files for many of the ISIS instruments can be found on the NMIDG webpage
- Follow the link [Moderator files for use with McStas](#)
- Download Commodus\_I.comp and then explore the Baseline and TS1 upgrade folders to see if the a suitable file exists. The TS1 instruments should all be there but only a few TS2 instruments have been added
- If not then the existing files distributed with McStas will work but will not have been updated to reflect the existing state of the moderators. If you want an updated file to be produced then you will need to ask the neutronics group to produce one.



# Using Commodus\_I

COMPONENT sourceMantid = Commodus\_I(  
Face="Larmor\_Base.mcstas", - data file in the folder with the .instr  
E0 = E\_min, - Minimum neutron energy  
E1 = E\_max, - Maximum neutron energy  
modXsize = 0.083, - moderator width  
modZsize = 0.033, - moderator height  
xw = 0.03, - target window width  
yh = 0.03, - target window height  
dist = 3.7) - distance to target window  
AT (0.0, 0.0, 0.0) RELATIVE Origin

Add this to the INITIALIZE section or E0=.....

E\_max=81.8047/(lmin\*lmin);

E\_min=81.8047/(lmax\*lmax);



# Concluding Comments

- Optimisation of an instrument is extremely complicated.
  - Define your metrics in discussion with others
  - Accept the fact that there will be a lot of compromises to be made.
  - Don't forget about misalignments.



# Exercises SANS

- SANS
  - Download the SANS2D instrument file from Github  
<https://github.com/McStasMcXtrace/McCode/tree/master/mcstas-comps/examples>
  - Compile and launch the trace so you can see what the instrument looks like.
  - Identify the key components of the instrument
    - Why is there a bender?
    - Why is the guide so small?
    - Why is the m value of the guide sections so low?
  - Modify the instrument to explore the incident divergence before the sample
  - Add a beamstop with a radius of 30mm to block the direct beam from hitting the detector
  - Modify the instrument to change the L1 and L2 and explore how this effects the q range
  - Modify the m value of the guides to see how this effects the incident flux
  - Add a dummy guide to explore the effect of gravity



# Exercises NR

- NR
  - Download the CRISP instrument file from Github  
<https://github.com/McStasMcXtrace/McCode/tree/master/mcstas-comps/examples>
  - Modify the file to use the new moderator component and compare the results with the version distributed with McStas. Take care to read the instructions detailing how to calculate the flux.
  - Modify the sample component to add layers to produce fringes.
  - Modify the instrument to allow the sample angle to be scanned.
  - Write a linear TOF detector or use monitor\_ND to write one.
  - Add the choppers and expore frame overlap effects

